

ÈLEC

PERFORM BETTER, COST LESS



SPEC

1:1
RATIO

CONTRA ANGLE SERIES CA100L(OPTIC) / CA100R(NON-OPTIC)



Model	CA100L
Ratio	1:1
Max.Speed	40,000RPM
Material	Stainless steel
Water	Internal water
Optic	Fiber Optic
Bur	CA bur(Ø2,35)



Model	CA100R
Ratio	1:1
Max.Speed	40,000RPM
Material	Stainless steel
Water	Internal water
Bur	CA bur(Ø2,35)

1:1
RATIO

STRAIGHT HANDPIECE SERIES SA100L(OPTIC) / SA100R(NON-OPTIC)



Model	SA100L
Ratio	1:1
Max.Speed	40,000RPM
Material	Stainless steel
Water	Internal water
Optic	Fiber Optic
Bur	HP bur(Ø2,35)



Model	SA100R
Ratio	1:1
Max.Speed	40,000RPM
Material	Stainless steel
Water	Internal water
Bur	HP bur(Ø2,35)

16:1
RATIO

ENDODONTIC HANDPIECE SERIES CA160L(OPTIC) / CA160(NON-OPTIC)



Model	CA160L
Ratio	16:1
Max.Speed	2,500RPM
Material	Stainless steel
Optic	Fiber Optic
Bur	CA bur(Ø2,35)



Model	CA160
Ratio	16:1
Max.Speed	2,500RPM
Material	Stainless steel
Bur	CA bur(Ø2,35)

20:1
RATIO

IMPLANT ANGLE SERIES SG200L(OPTIC) / SG200(NON-OPTIC)



Model	SG200L
Ratio	20:1
Max.Speed	2,000RPM
Max.Torque	70 N.cm
Material	Stainless steel
Optic	Fiber Optic
Bur	CA bur(Ø2,35)



Model	SG200
Ratio	20:1
Max.Speed	2,000RPM
Max.Torque	70 N.cm
Material	Stainless steel
Bur	CA bur(Ø2,35)